

TREATMENT FOR THE UPPER EXTREMITY

The Entré System, a segmented pneumatic device, provides therapy for patients with lymphedema by delivering intermittent compression to targeted extremities.

UPPER EXTREMITY TREATMENT SOLUTION



Control Unit (E0651)

The control unit is a non-calibrated, segmented pneumatic device which delivers compressed air via a hose connector that is attached to the extremity garment. The device applies gradient distal to proximal pressure. The control unit is programmable and can apply pressures ranging from 19 mmHg – 60 mmHg.

Upper Extremity Garment (E0668)

The upper extremity garment, made of easy-to-clean nylon material, has eight chambers that inflate in a distal to proximal motion to effectively move excess fluid from the swollen extremity. Garments are available in two different lengths, and are designed with zippers that allow a range of circumference sizes to be fitted.

Product Specifications


FDA Status: Cleared to market (Class II in May 2015)

Control Unit: Weight 4 pounds

Dimensions 6"W x 8"H x 11"D

Garment Fabric: Nylon (This product is not made with natural rubber latex)

TREATMENT PROTOCOL



The Entré System's programmability allows the clinician to choose from pressures ranging from 19 to 60 mmHg to address individual patient needs. Recommended protocol is one-hour therapy sessions once or twice daily per extremity.

Full Arm Treatment

Provides 60 minutes of therapy utilizing the full arm garment.