

FINALLY, AN EASIER WAY...
TO EFFECTIVELY MANAGE
YOUR CHRONIC SWELLING



AIRWEAR™
COMPRESSION SYSTEM

Tactile
MEDICAL®

Comfortable wrap-and-go air compression to ease swelling.

Introducing the Airwear™ wrap system:
an easy way to get effective, consistent swelling
relief throughout your day.

FEATURES AND BENEFITS

- No more struggling with stockings—just wrap on and pump to prescribed level of compression
- Consistently stays in place, providing you with support throughout the day
- Does not bind or pinch, allowing for a comfortable and effective fit
- Breathable, lightweight fabric fits discreetly under professional and casual clothing
- 3 compression levels to effectively control swelling from mild to severe days

MICHELE H.



It's probably the best thing I've ever had
for standing on my feet all day.


THOMAS K.



They've made so much of a difference to me as far as my
comfort level. Just wearing them, my legs felt good.

Ready to explore how Airwear can help you?

1. TAKE LEG MEASUREMENTS



- A** Mid-calf circumference measurement ____ cm
- B** Ankle circumference measurement ____ cm
- C** Mid-foot circumference measurement ____ cm

2. FIND YOUR SIZE AND COMPRESSION LEVEL

SIZE	A CALF	B ANKLE	C FOOT
Small	31–36 cm (≈12–14")	20–26 cm (≈8–10")	17–27 cm (≈6–10")
Medium	35–41 cm (≈14–16")	23–31 cm (≈9–12")	19–30 cm (≈8–12")
Large	40–48 cm (≈16–18")	28–37 cm (≈11–14")	22–33 cm (≈9–13")
X-Large	46–57 cm (≈18–22")	32–43 cm (≈13–17")	24–36 cm (≈10–14")

COMPRESSION LEVEL	20–30 mmHg	30–40 mmHg	40–50 mmHg
-------------------	------------	------------	------------

3. ORDER NOW

store.tactilemedical.com

or

833.924.7932 (833.9AIRWEAR)

IF YOU NEED AT-HOME MANAGEMENT
OF SWELLING ABOVE THE KNEE,
LEARN ABOUT THE CLINICALLY PROVEN
FLEXITOUCH® PLUS SYSTEM AT
TACTILEMEDICAL.COM

Tactile Medical is a leader in developing and marketing at-home therapy devices that treat lymphedema and chronic venous insufficiency.

Individual results may vary.

Indications/contraindications: Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

Tactile Medical

3701 Wayzata Blvd, Suite 300
Minneapolis, MN 55416 USA

Toll Free Tel: 833.924.7932 (833.9AIRWEAR)
Email: airwear@tactilemedical.com
Hours: 7 a.m. to 7 p.m. CT, Monday–Friday
store.tactilemedical.com

Tactile
MEDICAL®
HEALING RIGHT AT HOME

Flexitouch is a registered trademark and Airwear and Healing Right at Home are trademarks of Tactile Medical.
©2020 Tactile Medical. All rights reserved. 700679-001 Rev B 1/2020